


# KERN AET 200-4NM



Premium touchscreen analytical balance with the complete range of functions for demanding processes



## Functions

Capacity display:	yes
Conversion function:	yes
Counting function:	yes
Density determination function:	yes
Dispensing function:	yes
Function to save totals:	yes
Help function on the display:	yes
Percentage determination:	yes
PRE-TARE function:	yes
Statistical function:	yes
Time display:	yes
Vibration-free weighing (Animal weighing program):	yes
Weighing with tolerance range:	yes

## Environmental conditions

Maximum humidity:	80 %
Maximum operating temperature:	30 °C
Minimum ambient temperature:	18 °C

## Power supply

Input voltage:	110 V - 230 V AC
n/a:	EURO UK US

## Service

DAkkS Certificate:	963-101
Verification:	965-201

## Category

Brand:	KERN
Category:	Balances
Product Group:	Analytical balance

## Packaging & shipping

Delivery:	504 h
Dimensions packaging (WxDxH):	585 x 395 x 555 mm
Gross weight:	12,110 kg
Net weight:	7,134 kg
Shipping method:	Parcel service

## Construction

Cable length of display device:	35 cm
Casing material:	Plastic
Dimensions draft shield (WxDxH):	172 x 172 x 230 mm
Dimensions housing (WxDxH):	254 x 270 x 350 mm
Dimensions of display device (WxDxH):	215 x 156 x 71 mm

## Measuring system

Adjustment options:	Internal calibration (auto)
Linearity:	0,300 mg
n/a:	400 mg
Readout [d]:	0,100 mg
Reproducibility:	0,200 mg
Resolution:	2.200.000
Stabilisation time:	2 s
Warm up time:	8 h
Weighing range [Max]:	220 g
Weighing system:	Force compensation
Weighing units:	ct g mg

## Display

Display digit height:	2,100 cm
Display digit height:	1,900 cm

## Technical data about verification

Minimum load [Min]:	10 mg
Periodic Verification for Italy:	969-112
Verification approval conforming to 2014/31/EU:	yes
Verification by KERN possible:	yes
Verification class:	I
Verification value [e]:	1 mg

## Counting

Counting resolution:	2200000
Minimum piece weight at piece counting (Laboratory):	0,100 mg
Minimum piece weight at piece counting (Normal):	1 mg

# KERN AET200-4NM



Premium touchscreen analytical balance with the complete range of functions for demanding processes

Level indicator:	yes
Material weighing plate:	stainless steel
Overall dimensions mounted (WxDxH):	206 x 426 x 350 mm
Revolving screw feet:	yes
Weighing space (WxDxH):	168 x 160 x 227 mm
Weighing surface (d):	85